



Indication: For fixation around joints and trunks which are difficult to bandage. Can also be used for fixation of primary dressings and tubings. Only to be used on undamaged skin.

Fixation bandage

MEDIPLAST offers DANAMULL fixation bandages. The bandages are blood and salves repellent and have good aeration.

- Double weaving with woven edges.
- Elasticity: 100%
- Viscose, Polyamide
- Length: stretched

Art.No	Description	Ward box	Trsp.box
6050601	4 cm x 4 meters	20	1440
6050602	6 cm x 4 meters	20	960
6050603	8 cm x 4 meters	20	720
6050604	10 cm x 4 meters	20	600
6050605	12 cm x 4 meters	20	480

Manufactured in Denmark

Ward box = min order qty