


MEDI PLAST

Danafast DS

Tubular retention bandage, double stretch. Made of bamboo.



LONG-TERM FRESHNESS



WONDERFULLY SOFT



EXCELLENT HEAT REGULATION



ECO-FRIENDLY ✪

* <https://www.guaduabamboo.com/blog/environmental-impact-of-guadua-bamboo>
<https://www.trvst.world/sustainable-living/bamboo/all-about-bamboo-and-sustainability/>

MEDIPLAST® offers DANAFAST Double Stretch, a tubular elastic bandage that stretches in length and width and is thus easy to apply. Can be cut in desired length and can be applied in single or double layer, no additional fixation is required. Good air permeability. Colour coded with a woven line which makes choice of size easy.

TUBULAR BANDAGE DANAFAST DS BAMBOO

Bamboo. Length: 10 meters unstretched

Art.No	Description	Ward box	Trsp.box
6050785	RED = >3,5 cm, small arms and legs	1	35
6050786	GREEN = >5,5 cm, feet, hands, normal arms, small legs	1	28
6050787	BLUE = >8 cm, normal legs, large arms, thighs	1	20
6050788	YELLOW = >12 cm, small trunks, large thighs	1	14
6050789	BEIGE = >18 cm, adult trunks	1	7
6050790	PURPLE = >20 cm, large adult trunk	1	7

Bamboo. Length: 1 meter unstretched

6050791	RED = >3,5 cm, small arms and legs	72	576
6050792	GREEN = >5,5 cm, feet, hands, normal arms, small legs	72	576
6050793	BLUE = >8 cm, normal legs, large arms, thighs	72	432
6050794	YELLOW = >12 cm, small trunks, large thighs	36	288
6050795	BEIGE = >18 cm, adult trunks	36	144
6050796	PURPLE = >20 cm, large adult trunk	36	144

OEKO-TEX®
 CONFIDENCE IN TEXTILES
STANDARD 100



Indication: Can be used as a fixation dressing, under plaster and for dermatological use. Particularly suitable over joints, amputations and severe bandaged areas. Only to be used on undamaged skin.